Implant manual — Manuel d'implantation
Anwenderhandbuch — Manuale d'impianto
Manual de implante

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Implant manual

 SORIN | REPLY DR

Rate responsive dual-chamber pacemaker


****Information from the manufacturer****

The Dplus mode, Rate acceleration and AV delay shortening features indicated here in this manual are available only in REPLY devices marked with the following serial number: xxxZUxxx

SORIN Reply DR

Reminder

Battery depletion

Leads connection


TABLE OF CONTENTS

1. General description.....	6
2. Indications.....	6
3. Contraindications	7
3.1. Generalities	7
3.2. Potential complications	7
4. Warnings to patient.....	8
4.1. Risks related to medical environment.....	8
4.2. Storage.....	10
4.3. Characteristics of the sterile package	10
5. Implant procedure	10
5.1. Before opening the package	10
5.2. Preparing the pocket.....	11
5.3. Choosing the type of lead	11
5.4. Measurement of thresholds at implant.....	11
5.5. Lead connection.....	12
5.6. Functioning test at implant	13
6. Special modes	14
6.1. Nominal mode (Emergency mode)	14
6.2. Standby mode	14
6.3. Magnet mode	15
6.4. Response in the presence of interference	15
7. Functions and parameters	16
7.1. Pacing modes	16
♦ DDTA(R) Mode	16
♦ DDTV(R) Mode.....	17
♦ DDTAV(R) Mode	17
7.2. Rates.....	17
♦ Basic rate.....	17
♦ Maximum tracking rate	18
♦ Rate hysteresis.....	18
7.3. Pacing and sensing parameters	18
♦ Pulse amplitude	18

◆ Pulse width	19
◆ Sensitivity	19
◆ Pacing and sensing polarities.....	19
7.4. Refractory periods.....	20
◆ Absolute refractory periods	20
◆ Committed period	20
◆ Automatic AV delay	20
◆ AVD Paced/Sensed Offset.....	21
7.5. Rate response.....	21
7.6. Therapeutic improvements	22
◆ SafeR mode.....	22
◆ Dplus mode	22
◆ Rate acceleration.....	23
◆ AV delay shortening	23
◆ Rest rate	23
◆ Rate smoothing	24
7.7. Managing atrial arrhythmias (AA)	25
◆ Preventing AA: Overdrive	25
◆ Preventing AA: Post extrasystolic pause suppression	25
◆ Preventing AA: Acceleration on PAC	26
◆ Protecting the ventricles against AA: Mode switching	26
7.8. Optimizing operation	27
◆ Automatic detection of implantation	27
◆ Autosensing.....	28
◆ Ventricular Auto-threshold.....	29
◆ Anti-PMT protection.....	29

8. Patient follow-up30

8.1. Follow-up recommendations	30
8.2. Elective Replacement Indicator (ERI)	30
8.3. Explant	31
8.4. Programmer tests	32
8.5. Interrogating the memories (AIDA+)	32
8.6. Pacemaker identification.....	33

9. Physical characteristics	34
10. Electrical characteristics	34
10.1. Current drain at beginning of life	35
10.2. Correspondence between magnet rate and battery impedance	35
10.3. Battery	36
10.4. Longevity	36
11. Programmable parameters	37
12. Non-programmable parameters	41
13. Warranty	42
13.1. Article 1: Terms of warranty	42
13.2. Article 2: Terms of replacement	44

1. GENERAL DESCRIPTION

REPLY DR is a DDDR-type dual-chamber rate-responsive pacemaker.

It is equipped with a physiological sensor (minute ventilation) and an accelerometer to allow adaptation of pacing to suit the patient's activity.

REPLY DR provides a range of highly effective therapeutic functions:

- Prevention of atrial arrhythmias
- Protection of ventricles against atrial arrhythmias (Fallback Mode Switching)
- Preservation of atrioventricular conduction (SafeR mode and Dplus mode)
- Preservation of circadian cardiac rhythm (Rest rate)
- Prevention of vasovagal syncope (acceleration)
- Automatic adjustment of amplitude of ventricular pacing (Auto-threshold)
- Automatic adjustment of atrial and ventricular sensitivity (Autosensing)
- Protection against Pacemaker-Mediated Tachycardia (Anti-PMT)

2. INDICATIONS

Pacing indications are provided by the American College of Cardiology and the American Heart Association (see "ACC/AHA/HRS 2008 Guidelines for Device-Based Therapy of Cardiac Rhythm Abnormalities", J Am Coll Cardiol 2008 May 27;51(21):e1-62).

3. CONTRAINDICATIONS

3.1. GENERALITIES

There are no known contraindications to the use of cardiac pacemakers as a therapeutic method for controlling heart rate. Nevertheless, the patient's medical and physical condition, as well as age, should be taken into consideration when choosing the pacemaker and leads.

The benefits of pacing to pediatric subjects have not been evaluated. Adverse interactions may occur between the patient's spontaneous rate and pacemaker functions.

3.2. POTENTIAL COMPLICATIONS

Complications may arise with any implanted pacing system. They may be related to the device itself:

- Early battery depletion
- Component failure
- Circuit inhibition, reversion to standby mode, or other failures due to electromagnetic interference
- Pectoral stimulation

Complications may be related to the pacing leads:

- Poor connection to the pacemaker leads
- Lead displacement, cardiac wall perforation, or tissue reaction at the myocardium-electrode interface
- Insulation fracture
- Conductor fracture

Medical complications may also arise:

- Infection
- Fluid accumulation at the implant site

- Casing migration
- Skin erosion by the pacemaker with possible protrusion of the casing

4. WARNINGS TO PATIENT

The patient should be warned of the potential risks of pacemaker malfunction if he/she is exposed to external magnetic, electrical, or electromagnetic signals.

Electrical equipment: Household electrical appliances do not affect the functioning of the pacemaker, providing they are insulated to current standards. However, patients should avoid using induction ovens and cookers.

Antitheft gates: Since antitheft devices at the entrance to stores are not subject to any safety standards, it is advisable to spend as little time as possible in their vicinity.

Work environment: The patient's work environment may be an important source of interference. In that case, specific recommendations may be required.

Cellular phones: Radiofrequency signals can interfere with the functioning of REPLY DR if the handset is placed too close to the pacemaker. It is advisable to maintain a minimum distance of 15 cm (6 inches) between the cellular telephone and the implanted device, when the telephone is turned on.

4.1. RISKS RELATED TO MEDICAL ENVIRONMENT

Pacemaker operation should be carefully monitored prior to and after any medical treatment during which an electrical current from an external source passes through the patient's body.

Magnetic Resonance Imaging: MRI is strictly contraindicated in cardiac pacemaker patients.

Electrocautery or diathermy device: Diathermy and electrocautery equipment should not be used. If such devices must be used: **1.** Keep it as far as possible from the pacemaker. **2.** Set it at minimum intensity. **3.** Use it briefly.

External defibrillation: If the patient must be defibrillated, it is recommended that the paddles be placed at least 10 cm from the pacemaker, along a line perpendicular to the axis formed by the line connecting the leads to the pacemaker.

Internal defibrillator: Use of the pacemaker is contraindicated in implantable defibrillator patients.

Radiation therapy: Avoid exposure to ionizing radiation. Betatrons are contraindicated. If high doses of radiotherapy cannot be avoided, the pacemaker should be protected from direct exposure, and pacemaker operation should be monitored continuously. Resulting damage may be immediately undetectable. If irradiation of tissues close to the implantation site is necessary, it is recommended that the pacemaker be moved. As a safety measure, temporary pacemaker back-up should be immediately available.

Lithotripsy: The risks associated with lithotripsy are limited if the cardiac pacemaker is not implanted in the abdominal position. However, in order to avoid all risk of ventricular or atrial fibrillation, shocks should be administered synchronously with pacing.

Ultrasound: The pacemaker should not be exposed to therapeutic levels of ultrasound.

Transcutaneous Electrical Nerve Stimulation (TENS): TENS may interfere with pacemaker function. If necessary, the following measures may reduce interference: **1.** Place the TENS electrodes as close together as possible and as far as possible from the pulse generator and leads. **2.** Monitor cardiac activity during TENS use.

4.2. STORAGE

The pacemaker is packaged in a sterile pack contained within a cardboard storage package. It is recommended that the pacemaker be stored at a temperature ranging from 0 °C to 50 °C.

The cardboard package protects the pacemaker from impact during shipping.

Devices subjected to excessive impact, such as those that have been dropped on a hard floor, should not be implanted. Any device subjected to such an impact should be returned to your Sorin CRM representative for examination.

Devices **MUST NOT** be interrogated and programmed within the vicinity of other devices.

4.3. CHARACTERISTICS OF THE STERILE PACKAGE

The pacemaker and its ratchet screwdriver are sterilized with ethylene oxide and hermetically sealed in a transparent dual package in accordance with international standards.

5. IMPLANT PROCEDURE

5.1. BEFORE OPENING THE PACKAGE

Before opening the package:

- Check the "Use Before" date printed on the labels on the box and on the sterile package. Pacemakers not implanted before the "use before" date should be returned to your Sorin CRM representative.

- Interrogate the device: if a warning is displayed or if battery impedance is too high (> 2 kilohms), do not implant the device and contact your Sorin CRM representative.

Devices **MUST NOT** be interrogated and programmed within the vicinity of other devices.

Also check the integrity of the sterile package. The sterility of the contents is no longer guaranteed if the package has been pierced or altered. If the pacemaker is no longer sterile, it should be returned in its packaging to your Sorin CRM representative. Any re-sterilization of the unit is at the discretion of Sorin CRM.

5.2. PREPARING THE POCKET

Because of the orientation of the pacemaker connector, it is advisable to create the pocket in a left pectoral position.

In its final position, the pacemaker should be no more than 4 - 5 cm below the skin surface.

5.3. CHOOSING THE TYPE OF LEAD

For optimal use of pacemaker functions, the atrial lead should be bipolar. If a unipolar atrial lead is used, the Autosensing function will not be available.

Use of a high polarization ventricular lead interferes with normal operation of the Auto-threshold function.

Connections: The pacing/sensing connectors are IS-1 compatible. *Low profile* type leads are not compatible.

5.4. MEASUREMENT OF THRESHOLDS AT IMPLANT

Pacing and sensing thresholds should be measured at implant.

Pacing thresholds: Acute thresholds should be lower than 1 V for a 0.5 ms pulse width, both in the ventricle and in the atrium. Chronic thresholds should be less than 3 V for a 0.5 ms pulse width.

Sensing thresholds: For proper ventricular sensing, the amplitude of the R-wave should be greater than 5 mV. For proper atrial sensing, the amplitude of the P-wave should be greater than 2 mV.

5.5. LEAD CONNECTION


It is imperative that each lead be properly connected to the corresponding pacemaker connector (the position of each connector is indicated on the case).

Caution: Only both distal inserts must be tightened.

To connect each lead, proceed as follows:

1. Insert the ventricular lead pin all the way into the appropriate pacemaker port (marked V) of the pacemaker.
2. Insert the screwdriver in the centre of the prominent screwdriver slot contour on the distal insert, then into the pre-inserted screw. The ventricular screwdriver slot is positioned above the axis of the ventricular lead and the atrial screwdriver slot is positioned below the axis of the atrial lead.
3. Tighten the screw until the screwdriver clicks. Check the tightness by pulling the lead (the lead should remain in the pacemaker when pulling to check the tightness).
4. Insert the atrial lead pin all the way into the appropriate port (marked A) of the pacemaker.
5. Proceed as in steps 2 and 3.

Caution: 1. Do not tighten the pre-inserted screws when there is no lead connector inserted (this could damage the connector block). 2. Do not loosen the screws before inserting the lead connector (subsequent risk of being unable to reinsert the screw). 3. Use only the screwdriver provided with the pacemaker. Maintain the screwdriver's shaft perpendicular to the plane of the pacemaker (see figure below).

Warning: Ensure that the screwdriver's tip is fully inserted in the setscrew; otherwise the screwdriver can damage the setscrew and prevent connection with or disconnection from the lead.

To ensure full insertion, push the screwdriver's hex tip smoothly into the setscrew until it reaches the bottom of the hex chamber in the screw, which could be felt as a solid metallic contact. Do not implant the pacemaker if there is no feeling of solid metallic contact. Do not implant the pacemaker if the wrench does not click when attempting to tighten the setscrew on the lead pin.

5.6. FUNCTIONING TEST AT IMPLANT

After placing the device in its pocket, check the pacemaker functioning using an ECG.

6. SPECIAL MODES

6.1. NOMINAL MODE (EMERGENCY MODE)

A safety function may be rapidly restored by **pressing the following button**:

In nominal mode, the pacemaker operates with the parameters underlined in the table of programmable parameters.

6.2. STANDBY MODE

In the event of self-test failure, the pacemaker switches **automatically** to a safe mode of operation known as standby mode. In most cases, normal pacemaker operation may be re-established (if the pacemaker is at a stabilized temperature of approximately 37°C) by validating the reinitialization function suggested by the programmer. If this is not the case, contact your Sorin CRM representative.

In standby mode, the pacemaker operates with the following parameters: VVI, 70 min⁻¹, 5 V, 0.5 ms, sensitivity 2.2 mV, unipolar.

Programming is then impossible and test modes, including magnet mode, are ineffective.

6.3. MAGNET MODE

When a magnet is applied, the pacemaker terminates the cycle in progress and operates using the following parameters:

Mode	DOO (if mode programmed: Dxx or SafeR) VOO (if mode programmed: Vxx) AOO (if mode programmed: Axx)
Rate	Magnet rate
Pulse amplitude	5 V (or greater if the programmed amplitude value is greater)
Pulse width	0.5 ms (or greater if the programmed pulse width value is greater)
AV delay	Rest AV delay

On exiting magnet mode, the device operates:

- for 6 cycles by pacing at the magnet rate, with an AV delay of 95 ms (if applicable) and with the programmed amplitude,
- then for 2 cycles with the programmed parameters, in asynchronous mode.

6.4. RESPONSE IN THE PRESENCE OF INTERFERENCE

If the pacemaker senses electrical noise at a frequency above 20 Hz, it switches to an asynchronous mode at the basic rate. The programmed mode is restored as soon as the noise is no longer detected.

7. FUNCTIONS AND PARAMETERS

7.1. PACING MODES

The pacing mode determines pacemaker operation when pacing and sensing. The meanings of the codes (e.g., DDD, AOO) are provided in Appendix DD to Standard EN 45502-2-1. The SafeR mode is described in the section "Therapeutic improvements." The DDTA, DDTV, and DDTAV modes are special modes developed by Sorin CRM (see description below).

Contraindications: 1. AAI(R) pacing is contraindicated where AV conduction is impaired. 2. Asynchronous pacing is contraindicated when there is competition between the patient's intrinsic rhythm and the pacemaker. 3. DDD pacing is contraindicated in patients with chronic atrial fibrillation or flutter, or in those with slow retrograde conduction, which may induce Pacemaker- Mediated Tachycardia.

◆ DDTA(R) Mode

Dual-chamber DDD(R) pacing with triggered atrial pacing.

Both chambers are paced at the basic rate (DDTA) or at the sensor-driven rate (DDTAR) in the absence of intrinsic activity. A sensed P wave triggers atrial pacing and ventricular pacing after the programmed AV delay. A sensed R wave inhibits ventricular pacing if it occurs during the AV delay and inhibits atrial and ventricular pacing if it occurs outside the AV delay.

In this mode, atrial pacing may be triggered while maintaining atrioventricular synchronization.

This mode is not intended to be programmed permanently.

◆ DDTV(R) Mode

Dual-chamber DDD(R) pacing with triggered ventricular pacing.

Both chambers are paced at the basic rate (DDTV) or at the sensor-driven rate (DDTVR) in the absence of intrinsic activity. A sensed P wave inhibits atrial pacing and initiates ventricular pacing after the programmed AV delay. A sensed R wave triggers ventricular pacing and inhibits atrial pacing if it occurs outside the AV delay.

This mode allows triggered ventricular pacing while maintaining atrioventricular synchronization.

This mode is not intended to be programmed permanently.

◆ DDTAV(R) Mode

Dual-chamber DDD(R) pacing with triggered pacing in both chambers.

Both chambers are paced at the basic rate (DDTAV) or at the sensor-driven rate (DDTAVR) in the absence of intrinsic activity.

A sensed P wave triggers atrial pacing and ventricular pacing after the programmed AV delay.

A sensed R wave triggers ventricular pacing and inhibits atrial pacing if it occurs outside the AV delay.

This mode is not intended to be programmed permanently.

7.2. RATES

◆ Basic rate

The basic rate is the pacing rate in the absence of sensed events. If the patient's spontaneous rate falls below this limit, the pacemaker takes over.

Special case: If rest rate or hysteresis is programmed, the lower limit of the intrinsic rate may reach a value below the basic rate.

Programming requirement: If the rate-responsive pacing function is ON, the lower limit of the basic rate is 50 min⁻¹.

◆ **Maximum tracking rate**

The maximum tracking rate determines the upper limit of the pacing rate.

Special case: In triggered mode (single or dual chamber), the device may exceed the maximum tracking rate and pace to the rate limit value.

◆ **Rate hysteresis**

Rate hysteresis determines the ability of the patient's natural sinus rhythm to decrease below the escape rate. After sensing, the pacemaker operates with an escape rate reduced by the programmed percentage of hysteresis.

Programming requirements: 1. If the Rate Response function is programmed, the escape rate cannot be below 50 min⁻¹. 2. When the basic rate is less than or equal to 45 min⁻¹, rate hysteresis is programmed to 0% and is not accessible.

7.3. PACING AND SENSING PARAMETERS

◆ **Pulse amplitude**

The pulse amplitude determines the voltage applied to the heart during pacing.

Recommendation: Use the “as shipped” value immediately after implant to avoid loss of capture due to a temporary rise in the pacing threshold.

◆ Pulse width

The pulse width determines the duration for which the pulse amplitude is applied to the heart during pacing.

Programming requirement: When the Auto-threshold function is programmed to "Monitor" or "Auto", the programmable ventricular pulse width is less than or equal to 0.5 ms.

◆ Sensitivity

The sensitivity determines the minimum signal sensing threshold. For example, a sensitivity of 1 mV will allow sensing of only signals with an amplitude greater than 1 mV. In order to detect a low-amplitude cardiac signal, a higher sensitivity (lower value) must be programmed. However, if the pacemaker detects extracardiac signals, a lower sensitivity (higher value) must be programmed.

Recommendations: 1. A safety margin of 2 to 3 times the amplitude of the recorded signal should be maintained. 2. Programming a high atrial sensitivity (value less than or equal to 0.6 mV) should be reserved for P-waves of very low amplitude, since this increases sensitivity to external interference. These values may also be programmed for patients likely to present paroxysmal episodes of atrial fibrillation.

Programming requirements: 1. In unipolar mode, atrial sensitivity values below 0.4 mV are not available. 2. When the atrial (or ventricular) Autosensing is programmed to Auto, the atrial (or ventricular) sensitivity is no longer programmable.

◆ Pacing and sensing polarities

Although mechanically configured as a bipolar pulse generator, the pacemaker may be programmed either to unipolar or bipolar configuration. Pacing and sensing polarity configuration can be programmed independently for the atrial and ventricular channels.

The risk of pectoral stimulation is lower with bipolar pacing. With bipolar sensing, the pacemaker is less sensitive to myopotentials and to external electromagnetic noise.

Programming requirement: The pacemaker cannot be programmed to bipolar configuration if a bipolar lead is not connected (measurement of the impedance determines the type of lead present at the time of bipolar configuration is programmed).

7.4. REFRACTORY PERIODS

◆ Absolute refractory periods

Absolute atrial and ventricular refractory periods are initiated by any sensed or paced ventricular or atrial event. These periods are composed of a minimal Refractory Period and a triggerable Refractory Period. The duration of the refractory periods lengthens automatically as needed.

PVAB special case: The Post Ventricular Atrial Blanking (PVAB) is an absolute refractory period in which the minimal value is programmable.

◆ Committed period

In dual-chamber pacing, atrial pacing may be sensed by the ventricular sensing circuits, and may inhibit ventricular pacing. The pacemaker is protected against this phenomenon (called Crosstalk) by the committed period.

◆ Automatic AV delay

The AV delay is the programmable time interval between atrial sensing or pacing and ventricular pacing.

The rest AV delay is used at the programmed basic rate. The exercise AV delay is used at the programmed maximum rate. Between basic rate and maximum rate, the AV delay is calculated at each cycle by the pacemaker on the basis of a linear relationship between the AV delay and the atrial rate.

◆ AVD Paced/Sensed Offset

An AVD Paced/Sensed Offset is a period added to the AV delay after atrial pacing.

Programming requirement: In Dplus mode, the AV delay plus the AVD paced/sensed offset between atrial and ventricular pacing cannot exceed 350 ms. In the other pacing modes, the AV delay plus the AVD paced/sensed offset between atrial and ventricular pacing cannot exceed 300 ms.

7.5. RATE RESPONSE

The Rate Response function makes it possible to adjust the pacing rate based on the patient's physical activity.

Sensors used: The pacemaker is equipped with two sensors which prevent rate increases due to artefact through cross-checking of information.. **1. Minute ventilation (MV):** this is calculated based on the measurement of transthoracic impedance by standard (unipolar or bipolar) intracardiac leads. It is used to obtain a physiological response proportional to the degree of exertion during periods in which the patient is active. **2. An accelerometer (G),** which measures changes in the patient's anteroposterior acceleration: a rapid response may thus be obtained at the start of exertion and the end of exertion can be detected immediately.

Examples: **1. Automatic rate response** ("RR Auto" option). The pacing rate is constantly adapted to the patient's physical activity. **2. Manual rate response** ("RR Fixed" option). In that case, manual programming of the patients physical activity determines the pacing rate applied, based on the sensor signal. **3. Rate response during fallback mode switching** (SafeR/DDIR, DDD/DDIR and Dplus/DDIR modes). Rate response is inactive if the patient is in sinus rhythm. It becomes active when atrial arrhythmia is sensed.

Warnings: Use the Rate Response function cautiously in cases of:
1. Severe coronary insufficiency. 2. Severe aortic stenosis. 3. Myocardial function compromised by undue accelerations of the pacing rate.

Programming at implant: "Learn" and "RR Auto" should not be programmed prior to implant so as not to distort autocalibration.

Surgical procedure: For safety reasons, it is preferable to deprogram the Rate Response function before any surgical procedure on the pacemaker patient.

Programming requirement: 1. The SafeR/DDIR, DDD/DDIR and Dplus/DDIR modes (*Rate response during fallback mode switching*) are accessible only if Fallback mode switching is programmed.

Note: To protect the patient from prolonged pacing at high rates, rate responsive pacing is deactivated if the patient is paced at a rate higher than 100 min⁻¹ for a large number of cycles during any 24-hour period.

7.6. THERAPEUTIC IMPROVEMENTS

◆ SafeR mode

SafeR mode is intended to minimise unnecessary ventricular pacing.

The pacemaker functions in AAI mode, and temporarily switches to DDD mode upon the occurrence of AVB III, AVB II, AVB I and Pause.

◆ Dplus mode

Dplus pacing mode preserves spontaneous atrioventricular conduction.

Operation: An AV monitoring period allows small variations in the spontaneous atrioventricular delay. If no ventricular sensing occurs during that period, the pacemaker converts to DDD mode until atrioventricular conduction returns to normal. In DDD mode, the pacemaker uses the AV delays calculated from the averaged PR intervals.

The Dplus-R mode is a Dplus mode with Rate Response.

The Dplus/DDIR mode is a Dplus mode with Rate Response only during Mode Switching.

◆ **Rate acceleration**

The Acceleration function is used to suppress the cardiac pause and to reduce the vasodepressive response that accompanies it.

Indication: Acceleration should be used only for patients presenting sinocarotid hypersensitivity with vasodepressive response.

Operation: 1. The pacing rate is gradually increased (up to the programmed percentage) when an abrupt drop in sinus rate is sensed. 2. The rate is then gradually reduced in accordance with the principle of rate smoothing.

Limitation: The pacing rate induced by rate acceleration is limited to 120 min⁻¹ (or to the maximum rate if the latter is lower).

Programming constraint: 1. This function is available when Dplus or DDD mode is programmed. 2. If any of the AA prevention functions is activated, acceleration is forced to 0% and is not programmable.

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The Acceleration function (only) may be accompanied by a shortening of the AV delay, intended to increase hemodynamic response.

Indication: Shortening of the AV delay should be used only for patients presenting sinocarotid hypersensitivity with vasodepressive response.

◆ **Rest rate**

The rest rate is the pacing rate to which the pacemaker adjusts when it senses that the patient is sleeping or resting.

Sensing principle: Sleep or rest phases are determined by reduced respiratory and cardiac activity and a low occurrence of premature complexes.

Programming requirements: **1.** To be taken into consideration by the pacemaker, the rest rate must be below the basic rate. **2.** This function is programmable when the Rate Response mode is programmed to "RR Auto", "RR Fixed" or "Learn" and if sensor is "MV" or "MV + G".

◆ Rate smoothing

The Smoothing function is designed to prevent a sharp rate drop to the programmed basic rate in patients presenting with episodes of paroxysmal bradycardia.

Operation: **1.** If the patient's spontaneous rate drops, the pacemaker takes over at a rate slightly lower than the spontaneous rate immediately prior to the pause. **2.** The pacing rate is then gradually reduced until the basic rate is reached or the spontaneous rhythm is restored.

Programming at implant: It is not advisable to program the Smoothing function before implantation, since the pacemaker may detect noise and pace at a rate higher than the programmed basic rate.

Recommendation: Concomitant use of the Hysteresis function is highly advisable to avoid triggering of pacing for slight variations in the patient's spontaneous cycle (unless one of the AA prevention functions is activated).

Programming requirements: **1.** In VVI and VVT mode, the Smoothing function cannot be programmed to "Slow" or "Very slow" when the Auto-threshold function is programmed to "Monitor" or "Auto". **2.** The Smoothing function is automatically activated when an AA prevention function is activated.

7.7. MANAGING ATRIAL ARRHYTHMIAS (AA)

◆ Preventing AA: Overdrive

The Overdrive function is a function that prevents atrial arrhythmias. It is designed to pace the atrium as often as possible, while remaining close to sinus rhythm, within the limits of the programmed maximum Overdrive rate.

Contraindications: The use of this function is contraindicated if heart function is compromised by undue accelerations in the paced rate, and particularly in the event of: **1.** Severe coronary insufficiency. **2.** Severe aortic stenosis. **3.** Frequent ventricular arrhythmia. **4.** Advanced cardiomyopathy.

Side effect: Palpitations are possible. This function may actually induce high pacing rates after a premature atrial complex.

Programming constraint: When this function is programmed, the Fallback mode switching and Smoothing functions are automatically activated and the Acceleration and Hysteresis functions are automatically deactivated.

◆ Preventing AA: Post extrasystolic pause suppression

The Post extrasystolic pause suppression function is a function that prevents atrial arrhythmia.

Operation: This function is active in the event of a PAC or PVC. **1. Isolated PAC and short bursts.** In order to amortize the atrial interval following a PAC, an intermediate atrial escape interval is recycled at the PAC. **2. Suppression of the post PVC pause.** At a PVC, synchronous atrial pacing is triggered and a smoothed escape interval of 31 ms is recycled.

Warning: In case of myopotential or sensing of a T wave considered as a PVC, this function can induce an inappropriate acceleration.

Contraindications: The use of this function is contraindicated if heart function is compromised by undue accelerations in the paced rate, and particularly in the event of: **1.** Severe coronary insufficiency. **2.** Severe aortic stenosis. **3.** Frequent ventricular arrhythmia. **4.** Advanced cardiomyopathy.

Programming constraint: When this function is programmed, the Fallback mode switching and Smoothing functions are automatically activated and the Acceleration and Hysteresis functions are automatically deactivated.

◆ **Preventing AA: Acceleration on PAC**

The Acceleration on PAC function is a function that prevents atrial arrhythmia, intended to reduce occurrences of Premature Atrial Complexes (PAC). On each frequent isolated PAC, the pacing rate is increased by approximately 5 min^{-1} , within the rate limits calculated by the device.

Contraindications: The use of this function is contraindicated if heart function is compromised by undue accelerations in the paced rate, and particularly in the event of: **1.** Severe coronary insufficiency. **2.** Severe aortic stenosis. **3.** Sustained ventricular arrhythmia. **4.** Advanced cardiomyopathy.

Side effect: Palpitations are possible. This function may actually induce high pacing rates after a premature atrial complex.

Programming constraint: When this function is programmed, the Fallback mode switching and Smoothing functions are automatically activated and the Acceleration and Hysteresis functions are automatically deactivated.

◆ **Protecting the ventricles against AA: Mode switching**

The Mode switching function is designed to avoid prolonged ventricular pacing at a high rate for the entire duration of sustained atrial arrhythmia, by conversion to DDI(R) mode.

Operation: 1. Arrhythmia is identified by detection of abrupt acceleration of atrial rhythm. 2. When atrial arrhythmia is confirmed, the pacemaker converts to DDI(R) mode at a pacing rate that decreases gradually until it reaches the sensor-driven rate or the basic rate or the rest rate. 3. As soon as the arrhythmia ceases, the pacing rate is gradually increased up to the atrial rate (or sensor-driven rate). The pacemaker then converts to DDD mode: the atrium and ventricle are resynchronized.

Limitation: Mode switching occurs only if atrial arrhythmia is greater than 120 min⁻¹.

7.8. OPTIMIZING OPERATION

◆ Automatic detection of implantation

Operation: At implantation, once the ventricular lead has been connected to the pacemaker, 1. the device paces and senses the ventricle in both unipolar and bipolar fashion (Uni+Bi) by means of a temporary connection between the ring and the case. IF a bipolar ventricular lead is used AND the device is not in contact with the tissue, Uni+Bi pacing will result in bipolar pacing of the ventricle, which is a safety feature for pacing-dependant patients. 2. At the same time, the device automatically checks if the ventricular lead is still connected to the pacemaker and the case is in contact with the tissue by launching a unipolar ventricular impedance test. After five minutes, the device repeats this test to see if the ventricular lead is still connected to the pacemaker and the case still in contact with the tissue. If yes, the detection of implantation is confirmed. If no, a new check is done by the device every five minutes, until the ventricular lead connection and the contact between the pacemaker and the tissue are confirmed.

Polarity configuration: As soon as the detection of implantation has been confirmed, the device performs a polarity check by launching bipolar impedance tests for the atrial and ventricular leads. **1.** If the leads are confirmed to be bipolar, the sensing polarities are programmed to bipolar, and the pacing polarities are programmed to unipolar (as shipped values) or bipolar (if the values are re-programmed to "bipolar" in the box). **2.** If the leads are not confirmed to be bipolar, pacing and sensing are unipolar.

Measurements of lead impedances: Twenty minutes after the detection of implantation has been confirmed, the device starts measuring the lead impedances every 6 (six) hours and data is stored for each chamber.

Automatic launch of parameters: Twenty minutes after the detection of implantation: **1.** The device switches from DDD to SafeR mode. **2.** Memories are initiated. **3.** Sensor is programmed to "Learn". **4.** Statistics are reset to zero.

◆ Autosensing

The Autosensing function allows automatic adjustment of atrial and ventricular sensitivities.

Operation: **1.** If the atrial (or ventricular) Autosensing function is programmed to Auto, the atrial (or ventricular) sensitivity is constantly adjusted to 37.5% of the mean amplitude of the atrial (or ventricular) signals. **2.** When the atrium is paced, and for the duration of an atrial arrhythmia, the atrial sensitivity tends toward 0.4 mV. **3.** When the ventricle is paced, the ventricular sensitivity tends toward 1.5 mV in bipolar configuration and toward 2.5 mV in unipolar configuration.

Programming constraints: **1.** The atrial Autosensing function (Auto mode) is available only when the atrial sensing polarity is bipolar. **2.** When the Overdrive function is activated, the atrial Autosensing function is not programmable.

◆ Ventricular Auto-threshold

The Auto-threshold function allows automatic adjustment of the amplitude of ventricular pacing, according to a threshold test performed by the device at regular intervals.

Examples: **1.** If the Auto-threshold function is programmed to "Auto," the ventricular pacing amplitude is programmed to twice the pacing threshold. If the threshold search is not successful, the ventricular pacing amplitude is forced to 5 V. **2.** When the Auto-threshold function is programmed to "Monitor," the programmed ventricular pacing amplitude value is preserved.

Warning: Certain drugs can greatly increase the pacing threshold. The ventricular Auto-threshold function should be used with caution, properly adjusting the minimal ventricular amplitude.

Programming constraint: The ventricular Auto-threshold function is not programmable when the rest AV delay is less than or equal to 110 ms and/or when the basic rate is equal to 95 min⁻¹.

◆ Anti-PMT protection

Protection against Pacemaker-Mediated Tachycardia (PMT) can be provided for all patients with retrograde ventriculoatrial conduction without reducing atrial sensing capability.

Examples: **1. Palliative protection (permanently applicable).** The pacemaker does not trigger an AV delay on a P wave detected within 500 ms following an asynchronous ventricular event. **2. Termination mode ("TERMIN" option).** The pacemaker can identify these PMTs based on the stability of the VP interval (from ventricular paced beat to retrograde P wave). Once a PMT has been confirmed, the pacemaker lengthens the atrial refractory period in order to reduce the PMT. **3. Reprogramming mode ("REPROG" option).** In the event of recurrent PMTs, the pacemaker permanently shortens the rest and exercise AV delays.

Programming constraint: Automatic AV delay reduction is limited to 125 ms for rest AV delay and to 80 ms for exercise AV delay.

8. PATIENT FOLLOW-UP

8.1. FOLLOW-UP RECOMMENDATIONS

Standard follow-up: Annual standard follow up is recommended. Before the patient's discharge, and at each follow-up visit, it is advisable to: **1.** Check the battery status. **2.** Check for proper sensing (sensitivity, crosstalk) and pacing; set the pacing amplitude to twice the pacing threshold. **3.** Interrogate the device memories (AIDA+). **4.** Keep a printout of programmed parameters, test results, and memory data.

Frequent follow-up: It is advisable to follow up the patient every 6 months when the battery impedance becomes greater than or equal to 5 kilohms, especially for pacing-dependent patients.

Replacement: The pacemaker should be replaced as soon as the Elective Replacement Indicator (ERI) point is reached.

Implant software upgrade: In case a new implant software is downloaded in the device memory through the programmer, a warning message could be displayed by the programmer to inform the user and give the proper instructions to follow.

8.2. ELECTIVE REPLACEMENT INDICATOR (ERI)

The Elective Replacement Indicator (ERI) is controlled by:

- a magnet rate $80 \pm 1 \text{ min}^{-1}$
- battery impedance of 10 kilohms

Caution: The pacemaker should be replaced as soon as the Elective Replacement Indicator (ERI) point is reached.

When the ERI point is reached, the following parameters are switched to Nominal mode (listed in the table of programmable parameters) until the end of the device's life:

- Mode
- Basic rate
- Rate hysteresis
- Rate response mode
- Smoothing

Other parameters remain as programmed.

When the ERI point is reached and under following conditions: VVI, 70 min⁻¹, 0.4 ms, 500 Ohms, operation remaining before end of life is approximately:

- 9.9 months with 2.5 V programming
- 7.6 months with 3.5 V programming
- 5.9 months with 5 V programming
- 4.1 months with 7.5 V programming

The pacing amplitudes actually available at ERI are:

- 2.4 V with 2.5 V programming
- 3.45 V with 3.5 V programming
- 4.5 V with 5 V programming
- 6 V with 7.5 V programming

8.3. EXPLANT

The pacemaker should be explanted in the following cases:

- End of life
- Confirmed malfunction
- Cremation of the patient (the pacemaker may explode if placed in an incinerator)
- Patient death: for environmental reasons, the local regulation may require the explantation of the devices containing a battery supply.

The explanted pacemaker should not be reused in another patient.

All explanted pacemakers should be returned to your Sorin CRM representative carefully cleaned of all traces of contamination. Cleaning can be achieved by immersing the pacemaker in an aqueous sodium hypochlorite containing at least 1% chlorine, followed by rinsing copiously with water.

The pacemaker should be packaged so as to protect against the mechanical impact and temperature variations that may occur during shipping.

8.4. PROGRAMMER TESTS

The programmer is designed to:

- perform impedance tests on the atrial and ventricular leads,
- perform pacing threshold tests in asynchronous mode,
- determine the amplitude of sensed atrial and ventricular waves,
- display the patient's ECG with temporarily modified programming parameters (e.g., to observe the patient's subjacent rhythm),
- perform electrophysiologic studies (bursts and atrial and ventricular extrastimulus sequences).

8.5. INTERROGATING THE MEMORIES (AIDA+)

Interrogating the memories provides:


- the battery depletion curve,
- statistics since the last follow-up,
- 6-month curves and histograms for evaluating device operation,
- recorded episodes (fallback mode switching, atrial and ventricular bursts, switches from AAI to DDD) including intracardiac ECGs,
- diagnostic messages about events detected by the device.

8.6. PACEMAKER IDENTIFICATION

The pacemaker can be interrogated and programmed via telemetry, using the programming head interfaced with the Sorin CRM dedicated programmer.

The device can be non-invasively identified as follows:

1. Take an x-ray to identify the name of the manufacturer printed on the device (S = Sorin CRM ; P = Pacemaker ; K = Reply DR model).


2. Interrogate the device using the Sorin CRM dedicated programmer. The model and serial number of the device are automatically displayed. The first figure in the serial number corresponds to the last figure in the year of manufacture.

9. PHYSICAL CHARACTERISTICS

Dimensions	41.2 x 41.5 x 6.1 mm
Weight	20 g
Volume	8 cm ³
Connector	IS-1
Active surface area of casing	Non-coated model: 25.2 cm ² Coated model: 7 cm ²
Guaranteed connector holding power	14 N
Materials used	Active surface area of casing: 99% pure titanium Casing shell: silicone elastomer* Connectors: polyurethane*

* Medical-grade materials that have undergone "in vitro" and "in vivo" qualifications.

10. ELECTRICAL CHARACTERISTICS

Input impedance	Unipolar ≥ 40 kilohms Bipolar ≥ 40 kilohms
Pulse form	

10.1. CURRENT DRAIN AT BEGINNING OF LIFE

Shipment conditions (1)	safeR mode (3)	Inhibited: 6.5 μ A 100% pacing: 10.4 μ A
Shipment conditions (1)	DDD mode	Inhibited: 6.5 μ A 100% pacing: 14.2 μ A
Cenelec conditions (2)	safeR mode (3)	Inhibited: 6.9 μ A 100% pacing: 9.6 μ A
Cenelec conditions (2)	DDD mode	Inhibited: 6.9 μ A 100% pacing: 12.3 μ A

(1) 60 min⁻¹, 3.5 V, 0.35 ms, 500 Ohms, Holter ON, sensors ON.

(2) 70 min⁻¹, 2.5 V, 0.5 ms, 500 Ohms, Holter ON, sensors ON.

(3) Twenty minutes after the detection of implantation, the device switches automatically from DDD to safeR mode, assuming 1% ventricular pacing and 100% atrial pacing.

10.2. CORRESPONDENCE BETWEEN MAGNET RATE AND BATTERY IMPEDANCE

Magnet rate (min ⁻¹)	96	95	94	93	91	90
Magnet period (ms)	625	633	641	648	656	664
Battery resistance (kilohms)	<3	3	4	5	6	7
Magnet rate (min ⁻¹)	89	88	86	83	82	80
Magnet period (ms)	672	680	695	719	734	750
Battery resistance (kilohms)	7.5	8	8.5	9	9.5	10
Magnet rate (min ⁻¹)	78	75	73	70		
Magnet period (ms)	773	797	820	859		
Battery resistance (kilohms)	11	12	13	≥14		

10.3. BATTERY

Manufacturer	Greatbatch
Type	Lithium iodine
Model	GB8711
Total capacity	0.86 Ah
Usable capacity	BOL: 0.81 Ah. ERI: 0.05 Ah.
Voltage	BOL: 2.8 V. ERI: 2.5 V.

10.4. LONGEVITY

8.3 years	Shipment conditions (1) and safeR mode (3)
6.1 years	Shipment conditions (1) and DDD mode
9.0 years	Cenelec conditions (2) and safeR mode (3)
7.1 years	Cenelec conditions (2) and DDD mode

(1) 60 min^{-1} , 3.5 V, 0.35 ms, 500 Ohms, Holter ON, sensors ON.

(2) 70 min^{-1} , 2.5 V, 0.5 ms, 500 Ohms, Holter ON, sensors ON.

(3) Twenty minutes after the detection of implantation, the device switches automatically from DDD to safeR mode, assuming 1% ventricular pacing and 100% atrial pacing.

11. PROGRAMMABLE PARAMETERS

measured at 37 °C under a 500 ohm load

Legend:

Value in bold: "as shipped" value

Underlined value: nominal value

Basic parameters	Values
Mode ⁽¹⁾⁽²⁾	DDD -DDDR-DDD/DDIR-safeR- safeRR-safeR/DDIR-Dplus-Dplus- R-Dplus/DDIR-DDTA-DDTV- DDTAV-VDD-VDDR-DDI-DDIR- DOO- <u>VVI</u> -VVIR-VVT-VOO-AAI- AAIR-AAT-AOO-OOO
Basic rate (min-1) ⁽³⁾	30-40-45-50-55- 60 -65- <u>70</u> -75-80- 85-90-95 (± 2 min-1)
Rest rate (min-1)	50-55-60-65-70-75-80-85-90-95 (± 2 min-1)
Maximum tracking rate (min-1)	100-110-120- 130 -140-155-165- 175-185 (± 5 min-1)
Rate hysteresis (%)	<u>0</u> -5-10-20-35 (± 4 %)
Rest AV delay (ms)	30-45-65-80-95-110-125-140- 155 - 170-190-205-220-235-250 (± 10 ms)
Exercise AV delay (ms)	30-45-65- 80 -95-110-125-140-155- 170-190-205-220-235-250 (± 10 ms)
AVD Paced/Sensed Offset (ms)	0-15-30-45- 65 -80-95-110-125 (± 5 ms)

- (1) Twenty minutes after the automatic detection of implantation by the device, the as-shipped pacing mode (DDD) is automatically reprogrammed to SafeR.
- (2) Do not program OOO mode in pacing-dependent patients.
- (3) The corresponding basic periods and escape intervals are as follows:
1961-1500-1333-1200-1091-1000-923-857-800-750-706-667-632 ms.

Pacing/Sensing	Values
Atrial or ventricular amplitude (V) ⁽¹⁾	1.5-2.0-2.5-3.0- 3.5 -4.0- <u>5.0</u> -7.5 (± 20 %)
Atrial or ventricular pulse width (ms) ⁽¹⁾	0.10-0.25- 0.35 - <u>0.50</u> -0.60-0.75- 0.85-1 (± 35 µs)
Atrial sensitivity (mV) ⁽²⁾⁽³⁾⁽⁴⁾	0.1-0.2-0.3-0.4 (± 0.025 mV ± 40 %) 0.6-0.8- 1.0 -1.2-1.5-1.8-2.0-2.2-2.5- 2.7-3.0-3.5-4.0-4.5-5.0-6.0 (± 20 %)
Ventricular sensitivity (mV) ⁽²⁾⁽³⁾	1.0-1.2-1.5-1.8-2.0- <u>2.2</u> - 2.5 -2.7-3.0- 3.5-4.0-4.5-5.0-6.0-8.0-10.0-15.0 (± 20 %)
Atrial or ventricular sensing polarity ⁽⁵⁾	<u>Unipolar</u> -Bipolar
Atrial or ventricular pacing polarity ⁽⁵⁾	<u>Unipolar</u> -Bipolar

- (1) The nominal value may be greater than the underlined value (return to programmed value).
- (2) Values are measured using a positive and negative triangular signal of 2/13 ms.
- (3) For sensitivity settings below 2.5 mV, the pacemaker may detect noise lower than the level specified in clause 27.5.1 of standard EN 45502-2-1.
- (4) For sensitivity settings below 0.3 mV, the pacemaker may detect noise lower than the level specified in clause 27.4 of standard EN 45502-2-1
- (5) As soon as the detection of implantation has been confirmed, the lead configuration is automatically programmed to unipolar pacing and bipolar sensing (if a bipolar lead is used).

Special features	Values
Smoothing	OFF -Very slow-Slow-Medium-Fast
Acceleration (%)	0-5-15-25-35-45
AV delay shortening (ms)	0-15-30-45-65-80-95-110
Mode switching	ON -OFF
Anti-PMT protection	Termin- Reprog
Atrial or ventricular Autosensing	Auto- Monitor
Ventricular Auto-threshold	Auto-Monitor- OFF
Minimum ventricular amplitude (V)	1.5-2- 2.5 -3.0-3.5 ($\pm 10\%$)
Post V Atrial Blanking (ms) ⁽¹⁾	150 -165-180-195-210-225-240-255 (± 10 ms)

(1) After ventricular pacing, the minimal blanking in atrium is equal to the programmed value. After ventricular sensing, the minimal blanking in atrium is equal to the programmed value minus 55 ms.

AA prevention functions	Values
Pause suppression	A-V-A+V- OFF
Acceleration on PAC	ON -OFF
Overdrive	ON -OFF
Maximum Overdrive rate (min ⁻¹)	100-110-130-155-185 (± 15 min ⁻¹)

Rate response parameters	Values
Sensor choice	MV+G-MV-G
Rate response mode ⁽¹⁾	Learn-RR auto-RR fixed- <u>OFF</u>
Physical activity	Very low-Low-Medium-High-Very high
(1) Twenty minutes after the automatic detection of implantation by the device, the as-shipped rate response mode (OFF) is automatically reprogrammed to Learn.	
SafeR parameters	Values
Pause max. (s)	2- 3 -4
Long PR: max (ms)	250-300- 350 -400-450
Long PR: min (ms)	200- 250 -300-350-400-450
AVB I commutation	Rest+Exercise -Exercise
Automatic detection of implantation	Values
SafeR auto launch	Yes -No
Atrial pacing polarity	Unipolar -Bipolar
Ventricular pacing polarity	Unipolar -Bipolar

12. NON-PROGRAMMABLE PARAMETERS

Parameters	Values
Committed period (ms)	95 (± 5 ms)
Rate limitFast V Rate limit (min-1)	195 (± 5 min-1)
Refractory periods	automatic (see table below)

The minimum refractory period values are as follows (except in AAI, AAT, DDTA and DDTAV modes):

Event	Minimum atrial refractory periods (ms)	Minimum ventricular refractory periods (ms)
Atrium sensed	80 $\pm$ 10 ms	–
Atrium paced	AV delay	30 $\pm$ 2 ms
Ventricle sensed	programmable	95 $\pm$ 10 ms
Ventricle paced	programmable	150 $\pm$ 10 ms

The minimum refractory period values are in AAI, AAT, DDTA and DDTAV modes as follows:

Mode	Event	Minimum atrial refractory periods (ms)
AAI, AAT	Atrium sensed or paced	345 ms $\pm$ 10 ms
DDTA, DDTAV	Atrium sensed	205 ms $\pm$ 10 ms
DDTA, DDTAV	Atrium paced	AV delay

13. WARRANTY

The REPLY implantable cardiac pacemaker is the result of highly advanced research and all components have been selected after exhaustive testing.

Sorin CRM S.r.l. (identified as “Sorin CRM” hereafter) guarantees the product REPLY against any damage caused by component failure or production defects during a period of five years after the implantation date, and warranty commits itself to replace all REPLY devices according to the terms described in article 1 and described in article 2.

Sorin CRM makes no claim that the human body will not react unsuitably to the implantation of a REPLY, or that failure will never occur.

Sorin CRM does not guarantee the suitability of REPLY in defined types of patients: selection of the device is a medical decision.

Sorin CRM shall not be held liable for any damage indirectly associated with the REPLY, whether as part of normal or abnormal operation of the latter, nor damage from its explantation or replacement.

Sorin CRM does not authorise anyone to modify these warranty conditions.

13.1. ARTICLE 1: TERMS OF WARRANTY

1. The REPLY implantable cardiac pacemaker is only guaranteed for the first implantation.
2. The EURID/IAPM implant form must be sent to Sorin CRM within 30 days after implantation.
3. The REPLY cardiac pacemaker must be implanted prior to the use-before date indicated on the packaging.
4. The guarantee only applies to suspect devices returned to the manufacturer, carefully packed and accompanied by an explantation report duly completed by the hospital or the doctor and considered defective after analysis by Sorin CRM.

- The device must be returned within the 30 days following explantation to Sorin CRM.
 - Any device returned and replaced under the terms of this warranty will become the exclusive property of Sorin CRM.
 - Any rights under the terms of this warranty will be forfeited if the REPLY device has been opened by anyone other than Sorin CRM.
 - These rights will also be forfeited if the device has been damaged by carelessness or accident.
 - This is the case especially if the device has been exposed to temperatures above 50°C, to electrical abuse or to mechanical shock, particularly as a result of being dropped. Consequently, any expert opinion offered by a third party after the device has been removed also nullifies the guarantee.
5. The warranty will be forfeited if it is proven that the device has been misused or inadequately implanted, against the physicians' manual recommendations of REPLY.
 6. The warranty does not include leads and other accessories used for the implantation.
 7. The replacement terms or conditions described in article 2 include all devices that shall be replaced within the warranty period because of battery depletion, without any link to a component failure or a production hazard. The device battery longevity varies with the device settings and pacing operation over time.

13.2. ARTICLE 2: TERMS OF REPLACEMENT

1. In case of REPLY failure because of a component failure, a production defect, or a conception error, occurring within five-year period starting from the implantation date, Sorin CRM is committed to:
 - replacing free of charge the explanted device by a Sorin CRM device with equivalent features,
 - or issuing a replacement credit equal to the purchase price for the purchase of any other Sorin CRM replacement device.
2. In any case the credit issued by the warranty terms cannot exceed the purchase price of a Sorin CRM replacement device.

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