Implant manual – Manuel d'implantation
Anwenderhandbuch – Manuale di impianto
Manual de implante

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Implant manual

 SORIN | REPLY SR


Rate responsive single-chamber pacemaker

SORIN Reply SR

Reminder

Battery depletion

Lead connection


TABLE OF CONTENTS

1. General description.....	5
2. Indications.....	5
3. Contraindications	5
3.1. Generalities	5
3.2. Potential complications.....	6
4. Warnings to patient.....	6
4.1. Risks related to medical environment	7
4.2. Storage	8
4.3. Characteristics of the sterile package	9
5. Implant procedure	9
5.1. Before opening the package	9
5.2. Preparing the pocket	10
5.3. Choosing the type of lead.....	10
5.4. Measurement of thresholds at implant	10
5.5. Lead connection	10
5.6. Functioning test at implant	12
6. Special modes	12
6.1. Nominal mode (Emergency mode).....	12
6.2. Standby mode	12
6.3. Magnet mode.....	13
6.4. Response in the presence of interference.....	13
7. Functions and parameters	14
7.1. Pacing modes.....	14
7.2. Rates	14
♦ 7.2.1. Basic rate	14
♦ 7.2.2. Maximum tracking rate	14
♦ 7.2.3. Rate hysteresis	15
7.3. Pacing and sensing parameters.....	15
♦ 7.3.1. Pulse amplitude	15
♦ 7.3.2. Pulse width.....	15

♦ 7.3.3. Sensitivity	15
♦ 7.3.4. Pacing and sensing polarities	16
7.4. Refractory periods	16
♦ 7.4.1. Absolute refractory periods	16
7.5. Rate response	17
7.6. Therapeutic improvements	18
♦ 7.6.1. Rest rate	18
♦ 7.6.2. Rate smoothing	18
7.7. Optimizing operation	19
♦ 7.7.1. Automatic detection of implantation	19
♦ 7.7.2. Autosensing	20
♦ 7.7.3. Auto-threshold	20
8. Patient follow-up	21
8.1. Follow-up recommendations	21
8.2. Elective Replacement Indicator (ERI)	21
8.3. Explant.....	22
8.4. Programmer tests	23
8.5. Interrogating the memories (AIDA+)	23
8.6. Pacemaker identification	23
9. Physical characteristics	25
10. Electrical characteristics	25
10.1. Current drain at beginning of life	26
10.2. Correspondence between magnet rate and battery impedance	26
10.3. Battery	27
10.4. Longevity	27
11. Programmable parameters	28
12. Non-programmable parameters	30
13. Warranty	31
13.1. Article 1: Terms of warranty	31
13.2. Article 2: Terms of replacement	33

1. GENERAL DESCRIPTION

REPLY SR is a single-chamber rate-responsive pacemaker.

It is equipped with a physiological sensor (minute ventilation) and an accelerometer to allow adaptation of pacing to suit the patient's activity.

REPLY SR provides a range of highly effective therapeutic functions:

- Preservation of circadian cardiac rhythm (Rest rate)
- Automatic adjustment of amplitude of ventricular pacing (Auto-threshold)
- Automatic adjustment of sensitivity (Autosensing)

2. INDICATIONS

Pacing indications are provided by the American College of Cardiology and the American Heart Association (see "ACC/AHA/HRS 2008 Guidelines for Device-Based Therapy of Cardiac Rhythm Abnormalities", J Am Coll Cardiol 2008 May 27;51(21):e1-62).

3. CONTRAINDICATIONS

3.1. GENERALITIES

There are no known contraindications to the use of cardiac pacemakers as a therapeutic method for controlling heart rate. Nevertheless, the patient's medical and physical condition, as well as age, should be taken into consideration when choosing the pacemaker and leads.

The benefits of pacing to pediatric subjects have not been evaluated. Adverse interactions may occur between the patient's spontaneous rate and pacemaker functions.

3.2. POTENTIAL COMPLICATIONS

Complications may arise with any implanted pacing system. They may be related to the device itself:

- Early battery depletion
- Component failure
- Circuit inhibition, reversion to standby mode, or other failures due to electromagnetic interference
- Pectoral stimulation

Complications may be related to the pacing leads:

- Poor connection to the pacemaker leads
- Lead displacement, cardiac wall perforation, or tissue reaction at the myocardium-electrode interface
- Insulation fracture
- Conductor fracture

Medical complications may also arise:

- Infection
- Fluid accumulation at the implant site
- Casing migration
- Skin erosion by the pacemaker with possible protrusion of the casing

4. WARNINGS TO PATIENT

The patient should be warned of the potential risks of pacemaker malfunction if he/she is exposed to external magnetic, electrical, or electromagnetic signals.

Electrical equipment: Household electrical appliances do not affect the functioning of the pacemaker, providing they are insulated to current standards. However, patients should avoid using induction ovens and cookers.

Antitheft gates: Since antitheft devices at the entrance to stores are not subject to any safety standards, it is advisable to spend as little time as possible in their vicinity.

Work environment: The patient's work environment may be an important source of interference. In that case, specific recommendations may be required.

Cellular phones: Radiofrequency signals can interfere with the functioning of REPLY SR if the handset is placed too close to the pacemaker. It is advisable to maintain a minimum distance of 15 cm (6 inches) between the cellular telephone and the implanted device, when the telephone is turned on.

4.1. RISKS RELATED TO MEDICAL ENVIRONMENT

Pacemaker operation should be carefully monitored prior to and after any medical treatment during which an electrical current from an external source passes through the patient's body.

Magnetic Resonance Imaging: MRI is strictly contraindicated in cardiac pacemaker patients.

Electrocautery or diathermy device: Diathermy and electrocautery equipment should not be used. If such devices must be used:
1. Keep it as far as possible from the pacemaker. **2.** Set it at minimum intensity. **3.** Use it briefly.

External defibrillation: If the patient must be defibrillated, it is recommended that the paddles be placed at least 10 cm from the pacemaker, along a line perpendicular to the axis formed by the line connecting the leads to the pacemaker.

Internal defibrillator: Use of the pacemaker is contraindicated in implantable defibrillator patients.

Radiation therapy: Avoid exposure to ionizing radiation. Betatrons are contraindicated. If high doses of radiotherapy cannot be avoided, the pacemaker should be protected from direct exposure, and pacemaker operation should be monitored continuously. Resulting damage may be immediately undetectable. If irradiation of tissues close to the implantation site is necessary, it is recommended that the pacemaker be moved. As a safety measure, temporary pacemaker back-up should be immediately available.

Lithotripsy: The risks associated with lithotripsy are limited if the cardiac pacemaker is not implanted in the abdominal position. However, in order to avoid all risk of ventricular or atrial fibrillation, shocks should be administered synchronously with pacing.

Ultrasound: The pacemaker should not be exposed to therapeutic levels of ultrasound.

Transcutaneous Electrical Nerve Stimulation (TENS): TENS may interfere with pacemaker function. If necessary, the following measures may reduce interference: **1.** Place the TENS electrodes as close together as possible and as far as possible from the pulse generator and leads. **2.** Monitor cardiac activity during TENS use.

4.2. STORAGE

The pacemaker is packaged in a sterile pack contained within a cardboard storage package. It is recommended that the pacemaker be stored at a temperature ranging from 0 °C to 50 °C.

The cardboard package protects the pacemaker from impact during shipping.

Devices subjected to excessive impact, such as those that have been dropped on a hard floor, should not be implanted. Any device subjected to such an impact should be returned to your Sorin CRM representative for examination.

Devices **MUST NOT** be interrogated and programmed within the vicinity of other devices.

4.3. CHARACTERISTICS OF THE STERILE PACKAGE

The pacemaker and its ratchet screwdriver are sterilized with ethylene oxide and hermetically sealed in a transparent dual package in accordance with international standards.

5. IMPLANT PROCEDURE

5.1. BEFORE OPENING THE PACKAGE

Before opening the package:

- Check the "Use Before" date printed on the labels on the box and on the sterile package. Pacemakers not implanted before the "use before" date should be returned to your Sorin CRM representative.
- Interrogate the device: if a warning is displayed or if battery impedance is too high (> 2 kilohms), do not implant the device and contact your Sorin CRM representative.

Devices **MUST NOT** be interrogated and programmed within the vicinity of other devices.

Also check the integrity of the sterile package. The sterility of the contents is no longer guaranteed if the package has been pierced or altered. If the pacemaker is no longer sterile, it should be returned in its packaging to your Sorin CRM representative. Any re-sterilization of the unit is at the discretion of Sorin CRM.

5.2. PREPARING THE POCKET

Because of the orientation of the pacemaker connector, it is advisable to create the pocket in a left pectoral position.

In its final position, the pacemaker should be no more than 4 - 5 cm below the skin surface.

5.3. CHOOSING THE TYPE OF LEAD

For optimal use of pacemaker functions, the lead used should be bipolar. If a unipolar lead is used, the following will not be available:

- the Autosensing function in the atrial chamber,
- the rate-response function to minute ventilation.

Use of a high polarization lead interferes with normal operation of the Auto-threshold function.

Connections: The pacing/sensing connector is IS-1 compatible. *Low profile* type leads are not compatible.

5.4. MEASUREMENT OF THRESHOLDS AT IMPLANT

Pacing and sensing thresholds should be measured at implant.

Pacing thresholds: Acute thresholds should be lower than 1 V for a 0.5 ms pulse width. Chronic thresholds should be less than 3 V for a 0.5 ms pulse width.

Sensing thresholds: For proper ventricular sensing, the amplitude of the R-wave should be greater than 5 mV. For proper atrial sensing, the amplitude of the P-wave should be greater than 2 mV.

5.5. LEAD CONNECTION


It is imperative that the lead be properly connected to the pacemaker connector (the position of the connector is indicated on the case).

Caution: Only the distal insert must be tightened.

To connect the lead, proceed as follows:

1. Insert the lead pin all the way into the pacemaker port.
2. Insert the screwdriver in the centre of the prominent screwdriver slot contour on the distal insert, then into the pre-inserted screw.
3. Tighten the screw until the screwdriver clicks. Check the tightness by pulling the lead (the lead should remain in the pacemaker when pulling to check the tightness).

Caution: 1. Do not tighten the pre-inserted screw when there is no lead connector inserted (this could damage the connector). 2. Do not loosen the screw before inserting the lead connector (subsequent risk of being unable to reinsert the screw). 3. Use only the screwdriver provided with the pacemaker. Maintain the screwdriver's shaft perpendicular to the plane of the pacemaker (see figure below).

Warning: Ensure that the screwdriver's tip is fully inserted in the setscrew; otherwise the screwdriver can damage the setscrew and prevent connection with or disconnection from the lead.

To ensure full insertion, push the screwdriver's hex tip smoothly into the setscrew until it reaches the bottom of the hex chamber in the screw, which could be felt as a solid metallic contact. Do not implant the pacemaker if there is no feeling of solid metallic contact. Do not implant the pacemaker if the wrench does not click when attempting to tighten the setscrew on the lead pin.

5.6. FUNCTIONING TEST AT IMPLANT

After placing the device in its pocket, check the pacemaker functioning using an ECG.

6. SPECIAL MODES

6.1. NOMINAL MODE (EMERGENCY MODE)

A safety function may be rapidly restored by **pressing the following button**:

In nominal mode, the pacemaker operates with the parameters underlined in the table of programmable parameters.

6.2. STANDBY MODE

In the event of self-test failure, the pacemaker switches **automatically** to a safe mode of operation known as standby mode. In most cases, normal pacemaker operation may be re-established (if the pacemaker is at a stabilized temperature of approximately 37°C) by validating the reinitialization function suggested by the programmer. If this is not the case, contact your Sorin CRM representative.

In standby mode, the pacemaker operates with the following parameters: VVI or AAI (depending on the chamber), 70 min⁻¹, 5 V, 0.5 ms, sensitivity 2.2 mV, unipolar.

Programming is then impossible and test modes, including magnet mode, are ineffective.

6.3. MAGNET MODE

When a magnet is applied, the pacemaker terminates the cycle in progress and operates using the following parameters:

Mode	VOO (if mode programmed: Vxx) AOO (if mode programmed: Axx)
Rate	Magnet rate
Pulse amplitude	5 V (or greater if the programmed amplitude value is greater)
Pulse width	0.5 ms (or greater if the programmed pulse width value is greater)

On exiting magnet mode, the device operates:

- for 6 cycles by pacing at the magnet rate, with the programmed amplitude,
- then for 2 cycles with the programmed parameters, in asynchronous mode.

6.4. RESPONSE IN THE PRESENCE OF INTERFERENCE

If the pacemaker senses electrical noise at a frequency above 20 Hz, it switches to an asynchronous mode at the basic rate. The programmed mode is restored as soon as the noise is no longer detected.

7. FUNCTIONS AND PARAMETERS

7.1. PACING MODES

The pacing mode determines pacemaker operation when pacing and sensing. The meanings of the codes (e.g., VVI, AOO) are provided in Appendix DD to Standard EN 45502-2-1.

Contraindications: 1. AAI(R) pacing is contraindicated where AV conduction is impaired. 2. Asynchronous pacing is contraindicated when there is competition between the patient's intrinsic rhythm and the pacemaker.

7.2. RATES

◆ 7.2.1. Basic rate

The basic rate is the pacing rate in the absence of sensed events. If the patient's spontaneous rate falls below this limit, the pacemaker takes over.

Special case: If hysteresis is programmed, the lower limit of the intrinsic rate may reach a value below the basic rate.

Programming requirement: If the rate-responsive pacing function is ON, the lower limit of the basic rate is 50 min⁻¹.

◆ 7.2.2. Maximum tracking rate

The maximum tracking rate determines the upper limit of the pacing rate.

Special case: In triggered mode, the device may exceed the maximum tracking rate and pace to the rate limit value.

◆ 7.2.3. Rate hysteresis

Rate hysteresis determines the ability of the patient's natural sinus rhythm to decrease below the escape rate. After sensing, the pacemaker operates with an escape rate reduced by the programmed percentage of hysteresis.

Programming requirements: 1. If the Rate Response function is programmed, the escape rate cannot be below 50 min⁻¹. 2. When the basic rate is less than or equal to 45 min⁻¹, rate hysteresis is programmed to 0% and is not accessible.

7.3. PACING AND SENSING PARAMETERS

◆ 7.3.1. Pulse amplitude

The pulse amplitude determines the voltage applied to the heart during pacing.

Recommendation: Use the "as shipped" value immediately after implant to avoid loss of capture due to a temporary rise in the pacing threshold.

◆ 7.3.2. Pulse width

The pulse width determines the duration for which the pulse amplitude is applied to the heart during pacing.

Programming requirement: When the Auto-threshold function is programmed to "Monitor" or "Auto", the programmable ventricular pulse width is less than or equal to 0.5 ms.

◆ 7.3.3. Sensitivity

The sensitivity determines the minimum signal sensing threshold. For example, a sensitivity of 1 mV will allow sensing of only signals with an amplitude greater than 1 mV. In order to detect a low-amplitude cardiac signal, a higher sensitivity (lower value) must be

programmed. However, if the pacemaker detects extracardiac signals, a lower sensitivity (higher value) must be programmed.

Recommendations: 1. A safety margin of 2 to 3 times the amplitude of the recorded signal should be maintained. 2. Programming a high atrial sensitivity (value less than or equal to 0.6 mV) should be reserved for P-waves of very low amplitude, since this increases sensitivity to external interference.

Programming requirements: 1. When the Autosensing is programmed to Auto, the sensitivity is no longer programmable.

◆ 7.3.4. Pacing and sensing polarities

Although mechanically configured as a bipolar pulse generator, the pacemaker may be programmed either to unipolar or bipolar configuration. Pacing and sensing polarity configuration can be programmed independently.

The risk of pectoral stimulation is lower with bipolar pacing. With bipolar sensing, the pacemaker is less sensitive to myopotentials and to external electromagnetic noise.

Programming requirement: The pacemaker cannot be programmed to bipolar configuration if a bipolar lead is not connected (measurement of the impedance determines the type of lead present at the time of bipolar configuration is programmed).

7.4. REFRACTORY PERIODS

◆ 7.4.1. Absolute refractory periods

Absolute refractory periods are initiated by any sensed or paced event. These periods are composed of a minimal Refractory Period and a triggerable Refractory Period. The duration of the refractory periods lengthens automatically as needed.

7.5. RATE RESPONSE

The Rate Response function makes it possible to adjust the pacing rate based on the patient's physical activity.

Sensors used: The pacemaker is equipped with two sensors which prevent rate increases due to artefact through cross-checking of information. **1. Minute ventilation (MV):** this is calculated based on the measurement of transthoracic impedance by standard (bipolar) intracardiac lead. It is used to obtain a physiological response proportional to the degree of exertion during periods in which the patient is active.

2. An accelerometer (G), which measures changes in the patient's anteroposterior acceleration: a rapid response may thus be obtained at the start of exertion and the end of exertion can be detected immediately.

Examples: **1. Automatic rate response** ("RR Auto" option). The pacing rate is constantly adapted to the patient's physical activity. **2. Manual rate response** ("RR Fixed" option). In that case, manual programming of the patient's physical activity determines the pacing rate applied, based on the sensor signal.

Warnings: Use the Rate Response function cautiously in cases of: **1. Severe coronary insufficiency.** **2. Severe aortic stenosis.** **3. Myocardial function compromised by undue accelerations of the pacing rate.**

Programming at implant: "Learn" and "RR Auto" should not be programmed prior to implant so as not to distort autocalibration.

Surgical procedure: For safety reasons, it is preferable to deprogram the Rate Response function before any surgical procedure on the pacemaker patient.

Note: To protect the patient from prolonged pacing at high rates, rate responsive pacing is deactivated if the patient is paced at a rate higher than 100 min⁻¹ for a large number of cycles during any 24-hour period.

7.6. THERAPEUTIC IMPROVEMENTS

◆ 7.6.1. Rest rate

The rest rate is the pacing rate to which the pacemaker adjusts when it senses that the patient is sleeping or resting.

Sensing principle: Sleep or rest phases are determined by reduced respiratory and cardiac activity and a low occurrence of premature complexes.

Programming requirements: 1. To be taken into consideration by the pacemaker, the rest rate must be below the basic rate. 2. This function is programmable when the Rate Response mode is programmed to "RR Auto", "RR Fixed" or "Learn" and if sensor is "MV" or "MV + G".

◆ 7.6.2. Rate smoothing

The Smoothing function is designed to prevent a sharp rate drop to the programmed basic rate in patients presenting with episodes of paroxysmal bradycardia.

Operation: 1. If the patient's spontaneous rate drops, the pacemaker takes over at a rate slightly lower than the spontaneous rate immediately prior to the pause. 2. The pacing rate is then gradually reduced until the basic rate is reached or the spontaneous rhythm is restored.

Programming at implant: It is not advisable to program the Smoothing function before implantation, since the pacemaker may detect noise and pace at a rate higher than the programmed basic rate.

Recommendation: Concomitant use of the Hysteresis function is highly advisable to avoid triggering of pacing for slight variations in the patient's spontaneous cycle.

Programming requirement: 1. In VVI and VVT mode, the Smoothing function cannot be programmed to "Slow" or "Very slow" when the Auto-threshold function is programmed to "Monitor" or "Auto".

7.7. OPTIMIZING OPERATION

◆ 7.7.1. Automatic detection of implantation

Operation: At implantation, once the lead has been connected to the pacemaker, 1. the device paces and senses the chamber in both unipolar and bipolar fashion (Uni+Bi) by means of a temporary connection between the ring and the case. IF a bipolar lead is used AND the device is not in contact with the tissue, Uni+Bi pacing will result in bipolar pacing of the chamber, which is a safety feature for pacing-dependant patients. 2. At the same time, the device automatically checks if the lead is still connected to the pacemaker and the case is in contact with the tissue by launching a unipolar impedance test. After five minutes, the device repeats this test to see if the lead is still connected to the pacemaker and the case still in contact with the tissue. If yes, the detection of implantation is confirmed. If no, a new check is done by the device every five minutes, until the lead connection and the contact between the pacemaker and the tissue are confirmed.

Polarity configuration: As soon as the detection of implantation has been confirmed, the device performs a polarity check by launching a bipolar impedance test for the lead. 1. If the lead is confirmed to be bipolar, the sensing polarity is programmed to bipolar, and the pacing polarity is programmed to unipolar (as shipped values) or bipolar (if the value is re-programmed to "bipolar" in the box). 2. If the lead is not confirmed to be bipolar, pacing and sensing are unipolar.

Measurements of lead impedance: Twenty minutes after the detection of implantation has been confirmed, the device starts measuring the lead impedance every 6 (six) hours and data is stored.

Automatic launch of parameters: Twenty minutes after the detection of implantation: 1. Memories are initiated. 2. Sensor is programmed to "Learn". 3. Statistics are reset to zero.

◆ 7.7.2. Autosensing

The Autosensing function allows automatic adjustment of the sensitivity.

Operation: **1.** If the Autosensing function is programmed to Auto, the sensitivity is constantly adjusted to 37.5% of the mean amplitude of the sensed signals. **2.** When the rhythm is paced, the sensitivity tends toward 1.5 mV (ventricular chamber) or 0.4 mV (atrial chamber) in bipolar configuration and toward 2.5 mV (ventricular chamber) in unipolar configuration.

Programming constraint: **1.** The atrial Autosensing function (Auto mode) is available only when the atrial sensing polarity is bipolar.

◆ 7.7.3. Auto-threshold

The Auto-threshold function allows automatic adjustment of the amplitude of ventricular pacing, according to a threshold test performed by the device at regular intervals.

Examples: **1.** If the Auto-threshold function is programmed to "Auto," the ventricular pacing amplitude is programmed to twice the pacing threshold. If the threshold search is not successful, the ventricular pacing amplitude is forced to 5 V. **2.** When the Auto-threshold function is programmed to "Monitor," the programmed ventricular pacing amplitude value is preserved.

Warning: Certain drugs can greatly increase the pacing threshold. The ventricular Auto-threshold function should be used with caution, properly adjusting the minimal ventricular amplitude.

8. PATIENT FOLLOW-UP

8.1. FOLLOW-UP RECOMMENDATIONS

Standard follow-up: Annual standard follow up is recommended. Before the patient's discharge, and at each follow-up visit, it is advisable to: **1.** Check the battery status. **2.** Check for proper sensing (sensitivity) and pacing; set the pacing amplitude to twice the pacing threshold. **3.** Interrogate the device memories (AIDA+). **4.** Keep a printout of programmed parameters, test results, and memory data.

Frequent follow-up: It is advisable to follow up the patient every 6 months when the battery impedance becomes greater than or equal to 5 kilohms, especially for pacing-dependent patients.

Replacement: The pacemaker should be replaced as soon as the Elective Replacement Indicator (ERI) point is reached.

Implant software upgrade: In case a new implant software is downloaded in the device memory through the programmer, a warning message could be displayed by the programmer to inform the user and give the proper instructions to follow.

8.2. ELECTIVE REPLACEMENT INDICATOR (ERI)

The Elective Replacement Indicator (ERI) is controlled by:

- a magnet rate $80 \pm 1 \text{ min}^{-1}$
- battery impedance of 10 kilohms

Caution: The pacemaker should be replaced as soon as the Elective Replacement Indicator (ERI) point is reached.

When the ERI point is reached, the following parameters are switched to Nominal mode (listed in the table of programmable parameters) until the end of the device's life:

- Mode
- Basic rate

- Rate hysteresis
- Rate response mode
- Smoothing

Other parameters remain as programmed.

When the ERI point is reached and under following conditions: VVI (or AAI), 70 min⁻¹, 0.4 ms, 500 Ohms, operation remaining before end of life is approximately:

- 9.9 months with 2.5 V programming
- 7.6 months with 3.5 V programming
- 5.9 months with 5 V programming
- 4.1 months with 7.5 V programming

The pacing amplitudes actually available at ERI are:

- 2.4 V with 2.5 V programming
- 3.45 V with 3.5 V programming
- 4.5 V with 5 V programming
- 6 V with 7.5 V programming

8.3. EXPLANT

The pacemaker should be explanted in the following cases:

- End of life
- Confirmed malfunction
- Cremation of the patient (the pacemaker may explode if placed in an incinerator)
- Patient death: for environmental reasons, the local regulation may require the explantation of the devices containing a battery supply.

The explanted pacemaker should not be reused in another patient.

All explanted pacemakers should be returned to your Sorin CRM representative carefully cleaned of all traces of contamination. Cleaning can be achieved by immersing the pacemaker in an aqueous sodium hypochlorite containing at least 1% chlorine, followed by rinsing copiously with water.

The pacemaker should be packaged so as to protect against the mechanical impact and temperature variations that may occur during shipping.

8.4. PROGRAMMER TESTS

The programmer is designed to:

- perform impedance tests on the lead,
- perform pacing threshold tests in asynchronous mode,
- determine the amplitude of sensed waves,
- display the patient's ECG with temporarily modified programming parameters (e.g., to observe the patient's subjacent rhythm),
- perform electrophysiologic studies (bursts and extrastimulus sequences).

8.5. INTERROGATING THE MEMORIES (AIDA+)

Interrogating the memories provides:


- the battery depletion curve,
- statistics since the last follow-up,
- 6-month curves and histograms for evaluating device operation,
- diagnostic messages about events detected by the device.

8.6. PACEMAKER IDENTIFICATION

The pacemaker can be interrogated and programmed via telemetry, using the programming head interfaced with the Sorin CRM dedicated programmer.

The device can be non-invasively identified as follows:

1. Take an x-ray to identify the name of the manufacturer printed on the device (S = Sorin CRM ; P = Pacemaker ; N = Reply SR model).


2. Interrogate the device using the Sorin CRM dedicated programmer. The model and serial number of the device are automatically displayed. The first figure in the serial number corresponds to the last figure in the year of manufacture.

9. PHYSICAL CHARACTERISTICS

Dimensions	37.4 x 41.5 x 6.1 mm
Weight	19 g
Volume	7.5 cm ³
Connector	IS-1
Active surface area of casing	Non-coated model: 25.2 cm ² Coated model: 7 cm ²
Guaranteed connector holding power	14 N
Materials used	Active surface area of casing: 99% pure titanium Casing shell: silicone elastomer* Connectors: polyurethane*

* Medical-grade materials that have undergone "in vitro" and "in vivo" qualifications.

10. ELECTRICAL CHARACTERISTICS

Input impedance	Unipolar ≥ 40 kilohms Bipolar ≥ 40 kilohms
Pulse form	

10.1. CURRENT DRAIN AT BEGINNING OF LIFE

Shipment conditions (1)	VVI mode	Inhibited: 6.6 μA 100% pacing: 11.1 μA
-------------------------	----------	---

Cenelec conditions (2)	VVI mode	Inhibited: 6.6 μA 100% pacing: 9.3 μA
------------------------	----------	--

(1) 70 min^{-1} , 3.5 V, 0.35 ms, 500 Ohms, Holter ON, sensors ON.

(2) 70 min^{-1} , 2.5 V, 0.5 ms, 500 Ohms, Holter ON, sensors ON.

10.2. CORRESPONDENCE BETWEEN MAGNET RATE AND BATTERY IMPEDANCE

Magnet rate (min^{-1})	96	95	94	93	91	90
Magnet period (ms)	625	633	641	648	656	664
Battery resistance (kilohms)	<3	3	4	5	6	7
Magnet rate (min^{-1})	89	88	86	83	82	80
Magnet period (ms)	672	680	695	719	734	750
Battery resistance (kilohms)	7.5	8	8.5	9	9.5	10
Magnet rate (min^{-1})	78	75	73	70		
Magnet period (ms)	773	797	820	859		
Battery resistance (kilohms)	11	12	13	≥ 14		

10.3. BATTERY

Manufacturer	Greatbatch
Type	Lithium iodine
Model	GB8711
Total capacity	0.86 Ah
Usable capacity	BOL: 0.81 Ah. ERI: 0.05 Ah.
Voltage	BOL: 2.8 V. ERI: 2.5 V.

10.4. LONGEVITY

7.8 years	Shipment conditions (1) and VVI mode
9.3 years	Cenelec conditions (2) and VVI mode

(1) 70 min⁻¹, 3.5 V, 0.35 ms, 500 Ohms, Holter ON, sensors ON.

(2) 70 min⁻¹, 2.5 V, 0.5 ms, 500 Ohms, Holter ON, sensors ON.

11. PROGRAMMABLE PARAMETERS

measured at 37 °C under a 500 ohm load

Legend:

Value in bold: "as shipped" value

Underlined value: nominal value

Basic parameters	Values
Cavity	A-V
Mode ⁽¹⁾⁽²⁾	<u>VVI</u> -VVIR-VVT-VOO-AAI-AAIR-AAT-AOO-OOO
Basic rate (min-1) ⁽³⁾	30-40-45-50-55-60-65- <u>70</u> -75-80-85-90-95 (± 2 min-1)
Rest rate (min-1)	50-55-60-65-70-75-80-85-90-95 (± 2 min-1)
Maximum tracking rate (min-1)	100-110-120- 130 -140-155-165-175-185 (± 5 min-1)
Rate hysteresis (%)	<u>0</u> -5-10-20-35 (± 4 %)

- (1) The nominal value depends on the cavity selected.
- (2) Do not program OOO mode in pacing-dependent patients.
- (3) The corresponding basic periods and escape intervals are as follows:
1961-1500-1333-1200-1091-1000-923-857-800-750-706-667-632 ms.

Pacing/Sensing	Values
Pulse amplitude (V) ⁽¹⁾	1.5-2.0-2.5-3.0- 3.5 -4.0- <u>5.0</u> -7.5 (± 20 %)
Pulse width (ms) ⁽¹⁾	0.10-0.25- 0.35 - <u>0.50</u> -0.60-0.75-0.85-1 (± 35 µs)
Sensitivity (mV) ⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾	0.4 (± 0.025 mV ± 40 %) 0.6-0.8-1.0-1.2-1.5-1.8-2.0- <u>2.2</u> - 2.5 -2.7-3.0-3.5-4.0-4.5-5.0-6.0-8.0-10.0-15.0 (± 20 %)
Sensing polarity ⁽⁶⁾	<u>Unipolar</u> -Bipolar
Pacing polarity ⁽⁶⁾	<u>Unipolar</u> -Bipolar

- (1) The nominal value may be greater than the underlined value (return to programmed value).
- (2) Values are measured using a positive and negative triangular signal of 2/13 ms.
- (3) For sensitivity settings below 2.5 mV, the pacemaker may detect noise lower than the level specified in clause 27.5.1 of standard EN 45502-2-1.
- (4) Values below 1 mV are not available if the cavity selected is the ventricular cavity. Values above 6 mV are not available if the cavity selected is the atrial cavity.
- (5) If the cavity selected is the ventricular cavity, the “as-shipped” value is 2.5 mV. If the cavity selected is the atrial cavity, the “as-shipped” value is 1.0 mV.
- (6) As soon as the detection of implantation has been confirmed, the lead configuration is automatically programmed to unipolar pacing and bipolar sensing (if a bipolar lead is used).

Special features	Values
Smoothing	OFF -Very slow-Slow-Medium-Fast
Autosensing	Auto- Monitor
Ventricular Auto-threshold	Auto-Monitor- OFF
Minimum ventricular amplitude (V)	1.5-2- 2.5 -3.0-3.5 (± 10 %)

Rate response parameters	Values
Sensor choice	MV+G-MV-G
Rate response mode ⁽¹⁾	Learn-RR auto-RR fixed- OFF
Physical activity	Very low-Low-Medium-High-Very high
(1) Twenty minutes after the automatic detection of implantation by the device, the as-shipped rate response mode (OFF) is automatically reprogrammed to Learn.	
Automatic detection of implantation	Values
Pacing polarity	Unipolar-Bipolar

12. NON-PROGRAMMABLE PARAMETERS

Parameters	Values
Rate limitFast V Rate limit (min-1)	195 (± 5 min-1)
Refractory periods	automatic (see table below)

The minimum refractory period values are as follows:

Mode	Type of event	Minimum refractory periods (ms):
VVI, VVT	Ventricle sensed	95 ms (V chamber)
VVI, VVT	Ventricle paced	150 ms (V chamber)
AAI, AAT	Atrium sensed or paced	345 ms $\pm$ 10 ms (A chamber)

13. WARRANTY

The REPLY implantable cardiac pacemaker is the result of highly advanced research and all components have been selected after exhaustive testing.

Sorin CRM S.r.l. (identified as “Sorin CRM” hereafter) guarantees the product REPLY against any damage caused by component failure or production defects during a period of five years after the implantation date, and warranty commits itself to replace all REPLY devices according to the terms described in article 1 and described in article 2.

Sorin CRM makes no claim that the human body will not react unsuitably to the implantation of a REPLY, or that failure will never occur.

Sorin CRM does not guarantee the suitability of REPLY in defined types of patients: selection of the device is a medical decision.

Sorin CRM shall not be held liable for any damage indirectly associated with the REPLY, whether as part of normal or abnormal operation of the latter, nor damage from its explantation or replacement.

Sorin CRM does not authorise anyone to modify these warranty conditions.

13.1. ARTICLE 1: TERMS OF WARRANTY

1. The REPLY implantable cardiac pacemaker is only guaranteed for the first implantation.
2. The EURID/IAPM implant form must be sent to Sorin CRM within 30 days after implantation.
3. The REPLY cardiac pacemaker must be implanted prior to the use-before date indicated on the packaging.

4. The guarantee only applies to suspect devices returned to the manufacturer, carefully packed and accompanied by an explantation report duly completed by the hospital or the doctor and considered defective after analysis by Sorin CRM.
 - The device must be returned within the 30 days following explantation to Sorin CRM.
 - Any device returned and replaced under the terms of this warranty will become the exclusive property of Sorin CRM.
 - Any rights under the terms of this warranty will be forfeited if the REPLY device has been opened by anyone other than Sorin CRM.
 - These rights will also be forfeited if the device has been damaged by carelessness or accident.
 - This is the case especially if the device has been exposed to temperatures above 50°C, to electrical abuse or to mechanical shock, particularly as a result of being dropped. Consequently, any expert opinion offered by a third party after the device has been removed also nullifies the guarantee.
5. The warranty will be forfeited if it is proven that the device has been misused or inadequately implanted, against the physicians' manual recommendations of REPLY.
6. The warranty does not include leads and other accessories used for the implantation.
7. The replacement terms or conditions described in article 2 include all devices that shall be replaced within the warranty period because of battery depletion, without any link to a component failure or a production hazard. The device battery longevity varies with the device settings and pacing operation over time.

13.2. ARTICLE 2: TERMS OF REPLACEMENT

1. In case of REPLY failure because of a component failure, a production defect, or a conception error, occurring within five-year period starting from the implantation date, Sorin CRM is committed to:
 - replacing free of charge the explanted device by a Sorin CRM device with equivalent features,
 - or issuing a replacement credit equal to the purchase price for the purchase of any other Sorin CRM replacement device.
2. In any case the credit issued by the warranty terms cannot exceed the purchase price of a Sorin CRM replacement device.

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